















CLASS (OIML R111)	INDIVIDUAL STANDARD WEIGHTS	STANDARD WEIGHTS SETS		IMAGES
<b>E1</b>	1 mg – 50 kg	1 g – 50 g 1 g – 100 g 1 g – 200 g 1 g – 500 g 1 g – 1 kg 1 g – 5 kg 1 g – 10 kg	1mg – 50 g 1 mg – 100 g 1 mg – 200 g 1mg – 500 g 1 mg – 1 kg 1 mg – 2 kg 1 mg – 5 kg 1 mg – 10 kg	 
<b>E2</b>	1 mg – 50 kg	1 g – 50 g 1 g – 100 g 1 g – 200 g 1 g – 500 g 1 g – 1 kg 1 g – 2 kg 1 g – 5 kg 1 g – 10 kg	1mg – 50 g 1 mg – 100 g 1 mg – 200 g 1mg – 500 g 1 mg – 1 kg 1 mg – 2 kg 1 mg – 5 kg 1 mg – 10 kg	 
<b>F1</b>	1 mg – 50 kg	1 g – 50 g 1 g – 100 g 1 g – 200 g 1 g – 500 g 1 g – 1 kg 1 g – 2 kg 1 g – 5 kg 1 g – 10 kg	1mg – 50 g 1 mg – 100 g 1 mg – 200 g 1mg – 500 g 1 mg – 1 kg 1 mg – 2 kg 1 mg – 5 kg 1 mg – 10 kg	 
<b>F2</b>	1 mg – 50 kg	1 g – 50 g 1 g – 100 g 1 g – 200 g 1 g – 500 g 1 g – 1 kg 1 g – 2 kg 1 g – 5 kg 1 g – 10 kg	1mg – 50 g 1 mg – 100 g 1 mg – 200 g 1mg – 500 g 1 mg – 1 kg 1 mg – 2 kg 1 mg – 5 kg 1 mg – 10 kg	 
<b>M1</b>	1 mg – 2000 kg	1 g – 50 g 1 g – 100 g 1 g – 200 g 1 g – 500 g 1 g – 1 kg 1 g – 2 kg 1 g – 5 kg 1 g – 10 kg	1 mg – 50 g 1 mg – 100 g 1 mg – 200 g 1 mg – 500 g 1 mg – 1 kg 1 mg – 2 kg 1 mg – 5 kg 1 mg – 10 kg	 
<b>M2</b>	100 mg – 1000 kg	1 g – 50 g 1 g – 100 g 1 g – 200 g 1 g – 500 g 1 g – 1 kg 1 g – 2 kg 1 g – 5 kg 1 g – 10 kg	100 mg – 50 g 100 mg – 100 g 100 mg – 200 g 100mg – 500 g 100 mg – 1 kg 100 mg – 2 kg 100 mg – 5 kg 100 mg – 10 kg	 


<b>M3</b>	1 g – 2000 kg	1 g – 50 g	1g – 500 g	
		1 g – 100 g	1 g – 1 kg	
		1 g – 200 g	1g – 2 kg	
		1 g – 500 g	1g – 5 kg	
		1 g – 1 kg	1g – 10 kg	
		1 g – 2 kg		
1 g – 5 kg				
1 g – 10 kg				

## ACCESSORIES


TWEEZERS		
Type	Length (In mm)	Description
Size 1	105	Steel
Size 2	160	Steel
Size 3	230	Steel
Size 4	300	Steel
Size 1	105	Carbon with steel
Size 2	160	Carbon with steel
Size 3	230	Carbon with steel



GLOVES	
Type	Description
Size 1	Fine cotton gloves, thin quality
Size 2	Fine cotton gloves, strong quality
Size 3	Fine leather gloves with perforated cotton backhand



DUST BRUSHES	
Type	Diameter x high (in mm)
Size 1	D10 x L100
Size 2	D20 x L110
Size 3	D30 x L150
Size 4	D30 x L250



BOXES

Type

Plastic Boxes



Wooden Boxes

