



Sartorius ME235S



Incomparably fast. Incomparably stable accuracy. Incomparably reliable

- The new Sartorius ME takes only 8 seconds to deliver a stable and reproducible readout with an accuracy of 5 decimal places.
- Fully automatic draft shield system for fast adaptation to any work routine
- Built-in ionizer for neutralizing static electricity
- Ideal support when used as an inspection, measuring and test instrument in quality management systems
- isoCAL: fully automatic calibration/adjustment function with memory for calibration records
- SQmin function for monitoring minimum initial weight acc. to USP
- S.U.R.E. function for dynamic display of the uncertainty of measurement in accordance with DKD standards
- Three years warranty



Readability	0.01 mg
Weighing capacity	230 g
Repeatability	0.015 mg (60 g) 0.025 mg (60 to 230 g)
Max. linearity	0.1 mg
Response time (average)	8 s
Housing (WxDxH)	252 x 533 x 292 mm
Weighing pan	Ø 90 mm



Data printer



SartoConnect

Key accessories	Order No.
Data printer with date, time and statistics functions	YDP03-OCE
Foot switch for draft shield, tare and print functions	YPE01RC
Data transfer software SartoConnect with RS-232C standard cable (1.5 m);	YSC01L
Anti-static weighing pan	YWP01ME
Density determination kit	YDK01